

# V-100 Vision Screener



## Vision Screener V100


V100 is a portable diagnostic device to quickly detect vision abnormalities on patients

aged from 6 months to adults.

V100 can screen both eyes at once from 1 meter away.

Once the pupils are in focus, V100 will start to calculate automatically and deliver results

within 1 second.





### Quick Detection

- Myopia (nearsightedness)
- Anisometropia (unequal refractive power)
- Hyperopia (farsightedness)
- Strabismus (eye misalignment)
- Astigmatism (blurred vision)
- Anisocoria (unequal pupil size)

### Patient Management Software

MediView is the patient management system to help clinicians manage patient records and images and data collected by devices developed by MediWorks.

The testing result delivered by V100 can be quickly transferred to MediView via wifi connection.



Wi-Fi data transmission

Print via Bluetooth Printer



- Wireless data transfer by wifi
- Interesting sounds to attract children
- Super user-friendly with minimal user training required
- Capable of dealing with mass screening events

- Quickly deliver testing results in 1 second
- Automated screening provides thorough, objective, and easy-to-understand results
- Flip screen facilitate convenient vision measuring