


DERMATOME For Skin Transplantation




- M10 - 112**  
Power supply for direct corded use
- Mains voltage 100-240 Volts, 50-60 Hz
  - Supply voltage 15 Volts, current 4A
  - Protections: Short circuit / Overload / Over voltage
  - Dimensions: 140 x 128 x 89 mm
  - Weight: 560 g



- M10 - 113**  
Power cord for direct corded use
- Autoclavable,
  - Steam sterilisable at 2 bar 134° C
  - 4 m. cable length
  - 350 g



→ Cutting speed adjust button with safe lock.




- M10 - 110**  
Dermatome Handpiece
- Autoclavable
  - 16.000 cpm
  - 283 x 90 x 42 mm / 1080 g



→

Dermatome Cutting Thickness		
level	mm	inch
1	0.17	0.006
2	0.34	0.013
3	0.51	0.020
4	0.68	0.026
5	0.85	0.033
6	1.02	0.040
7	1.19	0.046




→ 8 flaps give the ability to adjust cutting width from 9.75 mm to 78 mm.



- M7 - 105**  
Spare blades  
One side cutting only.  
stainless steel
- 10 ea.
  - 81 x 18,5 x 0,40 mm

M10 - 111 Dermatome Electro Set	
M10-110	Dermatome Handpiece
M10-112	Power Supply
M10-113	Fixing socket with cable
M7-105	Spare blades x10 ea.
K6-247B	Sterile Container
M10-130	Basket with stand
K2-740	Paper filter x100 pcs
K2-660	Security seal x100 pcs
K2-672	Paper labels x100 pcs
K2-673	Paper labels x100 pcs



SET FOR STERILIZATION

- K2 - 740**  
Single use paper filter 100 pcs with indicator. Aygün mark.



- K2 - 673**  
Labels, paper x 100 pcs With indicator.



- M10 - 130**  
Basket with stand 405 x 253 x 70 mm



- K2 - 672**  
Labels, paper x 100 pcs With indicator.



- K6 - 247B**  
Sterile Container 481 x 280 x 135 mm



- K2 - 660**  
Security seal made of plastic 100 pcs

