

**insighters**  
M e d i c a l




**insighters Air-way Cowry**  
**We define the Innovation**

# Insighters Profile

Insighters Medical is an innovative High-Tech medical company located in Shenzhen, China. Its subsidiaries include 1 International Sales Center, 1 China Sales Center, 1 R&D Center, 2 Manufacturing Center, and 1 After-Sales Service & Quality Control Center.

Insighters has always been tightly understanding the demand from global clinical doctors, and strives to provide them the best airway management tools. The innovative Airway Cowry is a perfect solution for doctors in Anesthesia, ICU, Operation Room, Emergency Rescue to handle all kinds of airway management problems. With this, over 700 global hospitals become our customers, and fans of our innovation.

Insighters is in an open mind for various methods of cooperation with you. No matter you are a doctor, a dealer, a manufacturer, it's the same that you need a perfect airway management tool for patients, for hospitals, and for a mutual beneficial business, as we understand at its best quality, durability, and cost performance.



## insight iS3-C/F video flexible laryngoscope



OD	ID	ET tube	DLT	Oral Intubation	Nasal Intubation	Suction / Biopsy / Drug Delivery / Foreign Body remove / Oxygen supply
5.8mm	2.6mm	≥6.5	/	✓	/	✓
5.2mm	2.2mm	≥6.5	/	✓	✓	✓
4.5mm	1.2mm	≥5.5	F41	✓	✓	✓
3.6mm	/	≥4.5	F37	✓	✓	/
2.8mm	/	≥3.5	F32	✓	✓	/

By Insighters multi-function electronic flexible scope solution for OR/ICU/Respiratory department, doctors will face no troubles like hard operation, not visual or bad image view, worries of damaging the scopes easily like traditional Fiber optic bronchoscopes, etc.



▲  
**Wide Range  
of application**




▲  
**Fast intubation  
Clear and Reliable**

- Insight iS3-C/F is an innovative electronic flexible scope, without use any fiber-optic for imaging and illumination, it's not as easy to get damaged as traditional FOB, more durable for use and maintain.
- Insight iS3-C/F applies portable design, easy to move and carry for different departments use.
- Insight iS3-C/F has superb HD image quality, clear view for more accurate operation and diagnostics.



• Departments of Application •

Clear at a glance



Insighters Airway Management Set has a number of patentde technologies,it is the brilliant result of the team work and good communication with doctors, simple and powerful, which helps you to cope with the airway management with better understanding.

**Leading Air-way Management Visualization**



**Shenzhen Insighters Medical Technology Co., Ltd.**



Add: Floor 13, HengTengMei building, NO.3 Gan Li Road, Long Gang district, ShenZhen city, GaungDong Province, 518112,China.  
Tel: +86-755-28218609,28218551 Fax:+86-755-84509803 E-mail: [contact@insighters.cn](mailto:contact@insighters.cn) Web: [www.insighters.cn](http://www.insighters.cn)