

Balance ou moniteur pour poches à sang

Réf. SOC-BBE-001



- Balance ou moniteur pour poches à sang
- Blood Collection Monitor is a compact instrument to provide smooth and gentle rocking for homogeneous mixing with anticoagulant without clot formation of blood cells during collection of blood from a donor.
- Salient Features
 - Volume can be set in 1 ml. increments.
 - Provision of pausing collection and change.
 - Programmed volume during pause.
 - Micro-controller based program.
 - Volume can be set from 1 ml – 600 ml.
 - Display of weight and volume.
 - Auto tare facility.
 - Motor activated clamping at the end of the collection.
 - Audio - Visual alarm.
 - Auto Calibration.
 - Time totaliser.



Le fabricant et nous même, nous réservons le droit de modifier les produits en vue d'une amélioration