

electrophoresis automated densitometer

MiniSCAN



Pocket Electrophoresis Densitometer

- Ready for cleared/uncleared scanning
- Fast and reliable
- Suitable for agarose gel and cellulose acetate


socimed

VENTE DE CONSOMMABLES ET D'APPAREILS MÉDICAUX
Tel. +331 49 98 86 72 - socimed@socimed.com



ELECTROPHORESIS AUTOMATED DENSITOMETER

smart and reliable



allows to perform the most spread techniques: Serum Proteins, Haemoglobins, Lipoproteins

Technical features

- LED light source and high gain fototransistor
- Reading system in absorbance with self-regulation of the source bright
- Reading media in cellulose acetate cleared and uncleared
- Reading media in cellulose acetate and agarose gel
- 8 reading samples
- Connection to PC via USB
- External power Supply 24V
- Dimensions: 116x156x63 mm (LPH)

Software Features

- Pre-programmed Analysis: Serum proteins, Haemoglobins, Lipoproteins
- Possibility to define additional exams
- Automatic Processing paths electrophoretic pattern
- Possibility of correction of paths
- Personal Management and printing reports
- Report header customization
- Operating system: Microsoft Windows XP, Windows Vista and Windows 7

CE

IVD

...your partner in Protein Electrophoresis

socimed

VENTE DE CONSOMMABLES ET D'APPAREILS MÉDICAUX
Tel. +331 49 98 86 72 - socimed@socimed.com

