

support you care you



WELL SOLD FOR 40 YEARS WORLDWIDE
MEET THE BASIC DEMANDS OF GENERAL SURGERY

MT300 Manual Hydraulic Operating Table

MT300



Introduction

MT300 is widely used in chest, Abdominal surgery , ENT , gynaecology and obstetrics , urology and orthopedics etc.

Hydraulic lift by foot pedal, head operated movements.

The base and column cover all made by premium 304 stainless steel.

Tabletop is made of Composite laminate for X-ray, makes a High definition image.

support you care you

MT300 Manual Hydraulic Operating Table

WELL SOLD FOR 40 YEARS WORLDWIDE
MEET THE BASIC DEMANDS OF GENERAL SURGERY



Main Technical Data

Parameters may have reasonable error.

Item	MT300
Length and width for tabletop	2010/500mm
Elevation(down/up)	760/1010mm
Trendelenburg/anti-trendelenburg	22°/22°
Lateral tilt(left/right)	20°/20°
Head plate(up/down)	60°/70°
Leg Plate(up/down/outward)	15°/85°/90°
Back plate(up/down)	75°/15°
Kidney bridge up (Optional)	120mm