

U60

Diagnostic Ultrasound System

Clear. Compact. Comprehensive



EDAN

A world of potential

U60 Diagnostic Ultrasound System

With its compact design, user-friendly operation and advanced imaging technologies, the U60 diagnostic ultrasound system offers a broad range of clinical applications to address all of your imaging needs.

Compact Design

- Compact, lightweight design facilitates mobility
- Integrated battery enables operation anywhere
- Large 15" monitor for easy viewing

Comprehensive Clinical Applications


- Full set of transducers to support all applications

Advanced Imaging Technologies

- Tissue Harmonic Imaging
- Spatial Compounding Imaging
- Speckle Reduction Imaging


User-Friendly Operation

- Backlit, easy-to-use control panel
- Customizable user defined controls
- Large capacity data storage




Excellent Image Quality


The U60 provides excellent detail, contrast resolution and image uniformity without sacrificing penetration. Excellent color and Doppler sensitivity easily detects small, low flow vessels.




Common Carotid Artery




Right Hepatic Vein




Mitral Valve



Renal Vascularity



Interscalene



Umbilical Artery

Transducers



2-6MHz
Curved Array



2-6MHz
Curved Array

5-8MHz
Endocavity

2-4MHz Phased Array



7-16MHz
Linear Array



8-12MHz
Linear Array



5-10MHz
Linear Array



4-8 MHz
Linear Array



3-6MHz
Micro-Convex

4-9MHz
Micro-Convex

5-10MHz
Micro-Convex

About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- In-Vitro Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



Edan Instruments, Inc. | No.15 Jinhui Rd., Jinsha Community,
Kengzi Subdistrict, Pingshan District, Shenzhen | 518122, P.R. China

© 2017 Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all territories, contact Edan for local availability.



ENG-US-U60-V1.1-20180305